




PETIETEC

**PRODUCT INTRODUCTION
MANUAL FOR PET
HEATING PHYSIOTHERAPY
BRACES**

HIGH-TECH GRAPHENE HEATING



Upgraded Waterproof
3 in 1 Thermal Fabric



Easy to use
power bank

PRODUCT FEATURES AND PRINCIPLE INTRODUCTION

Features:


- Veterinarian Approved
- Targeted Heating & 3 Temperature options
- Made of Graphene material without over heating
- Totally safe with no sided effects
- Work well with other therapies

Principle of action:

- Hyperthermia/infrared rays relieve pain and speed up local blood circulation
- Heating improves the lubricating performance of joint lubricating fluid

WHO CAN BENEFIT FROM OUR HEATING THERAPY BRACES

- * Dogs with chronic pain, like arthritis and hip/elbow dysplasia
- * Dogs with Muscle Stiffness/Soreness
- * Dogs with a lot of exercise and training, like working dogs and agility dogs
- * Senior Dogs



OCCASIONS TO USE HEAT THERAPY



- 1.FOR RELIEVING JOINT OR MUSCLE PAIN AND STIFFNESS
- 2.MANAGE CHRONIC ARTHRITIS SYMPTOMS
- 3.RECOVERY FROM JOINTS SURGERY
- 4.WARM UP BEFORE EXERCISE OR REHABILITATION
- 5.SOOTH SORE MUSCLES AFTER A WORKOUT
6. TO TREAT MUSCLE INJURIES

TYPES OF BRACE



INSTRUCTIONS

1. Your dog may need a few minutes to adjust to wearing the braces. Check the fit to ensure there are no areas of discomfort. Monitor the dog during the treatment period.
2. Vet suggests doing the treatment 15-30 minutes each time, 2 times a day.
3. This product has 3 levels of temperature, and the highest temperature has been rigorously tested and verified. In addition, we use graphene heating material, and the heat conduction is very uniform and comfortable, safe and no any side effects.

The three temperature ranges:

Red – (48°C / 118.4°F)

Blue – (40°C / 104°F)

Green – (35°C / 95°F)

Scan QR code for manual video!



Heated therapy vest



Hip heated therapy brace



Knee/elbow heated therapy brace



Shoulder & Neck heated therapy brace

CUSTOMERS REVIEWS

1

Christelle Coste

FRANCE • Police dog trainer



I am a dog trainer in the municipal police police officer. Fenvers, my dog, an 11 year old dutch shepherd, retired police dog. He was a competition dog at Mondioring and he made two participations in the world championship. He is currently having difficulty to move properly and he suffers from joint pain in his pelvis and his back legs.

My dog tolerates the PetieTec orthosis very well :) After using it, he can go out again to exercise. To go out and do some exercises. I want to say that PetieTec is great I Thank you very much for creating a product a product that really works!

2

Emily

US • Veterinary physiotherapy





I am currently using them for a client of mine and wow!

I used the braces for a lab with severe joint issues from my client and my own cocker spaniel as working dog with quite tight lower back. I really am supportive of these products, I think they are really cool! I also use 10mins prior to exercise to warm up the muscles and 10 mins in the morning/after long periods of rest to help with any pain. After using it, he can go out again to exercise. To go out and do some exercises. I want to say that PetieTec is great I Thank you very much for creating a product a product that really works!

3

Aurélie Hrs

FRANCE • Animal Physiotherapist



Petietec Heating braces are very interesting for the follow-up of senior dogs or dogs with locomotion difficulties. In particular for dogs with difficulties to stand up cold.

The contribution of heat allows a better mobility as of the rise or daily during the day. The heat relieves pain. As an animal welfare professional, I use this product with pleasure during my sessions.